



evomed evoscope

Family of Endoscopic Products


EVOMED, offers a growing range of laparoscopic access devices that aim to meet the majority of diverse laparoscopic surgical needs. In the development of these products primary focus is on absolute safety to the patient achieved through ergonomic design and ease of use.

Please contact us for any additional information on any of the evomed products.

evoreach Endoscopic Instruments

The EVOREACH range of single patient use endoscopic instruments are intended to provide consistent and precise performance that assures the surgeon of safety and reliability. All instruments are 5mm in diameter, have rotating 33cm insulated shafts, all include a mono polar electrocautery feature that meets most of the needs of the endoscopic surgeon. All graspers have a lockable ratchet with quick release feature.

Design focus is on the ergonomic design of the handles to minimise operator hand fatigue and the working actuator jaws that are designed based on traditional reusable instrument technology, with maximum function as required, that is, precision jaw alignment and cutting (sharp to the blade tip), dissecting and grasping.



FEATURES	ADVANTAGES	BENEFITS
Ergonomic handle	Easier surgeon handling and surgeon comfort	Reduces hand fatigue
Wide jaw opening of 60°	Larger tissue bites Improved manipulation of tissue	More efficient tissue handling and procedure
Fully insulated shaft	Eliminates stray electrical current	Improves safety to surgeon and patients
360° Finger wheel rotation	Improves access to target tissue and precise jaw alignment	Ease of handling and reduce hand fatigue
Angled electosurgery terminal post	Removes cables away from surgical field	Provides unobstructed surgical field and organisation
Ratchet trigger on graspers	Easy locking and releasing of tissue grasping	Can engage or disengage ratchet
Uses mono polar electricity	Coagulation option through every instrument	Additional operating efficiency

PRODUCT IMAGES	PRODUCT REFERENCE	DESCRIPTION	DIAMETER (mm)	360° ROTATION	LENGTH (cm)	RATCHET
	AG-05	Atraumatic Fenestrated Grasper (Short), lockable with ratchet	5 mm	YES	33 cm	YES
	GA-05	Allis Grasping Forcep lockable with ratchet	5 mm	YES	33 cm	YES
	EW-05	Atraumatic Fenestrated Intestinal Grasper (Long), lockable with ratchet	5 mm	YES	33 cm	YES
	EB-05	Babcock Grasping Forcep, lockable with ratchet	5 mm	YES	33 cm	YES
	GJ-05	Standard Straight Grasper, lockable with ratchet	5 mm	YES	33 cm	YES
	EC-05	Endo Clinch Grasper 5mm, lockable with ratchet	5 mm	YES	33 cm	YES
	ED-05	Maryland Blunt Tip Dissecting Forcep	5 mm	YES	33 cm	NO
	ESC-05	Curved (Metzenbaum) Scissor	5 mm	YES	33 cm	NO



EVOMED Group Sarl

pH8, pépinière d'entreprises innovantes
8 place de l'Hôpital
67000 STRASBOURG, France
Phone: + 33 (0) 3 88 23 70 40
Fax: + 33 (0) 3 88 23 70 00
Web: www.evomedgroup.com



Please contact us for any additional information on any of the evomed products.