



evomed evoscope

Family of Endoscopic Products


EVOMED, offers a growing range of laparoscopic access devices that aim to meet the majority of diverse laparoscopic surgical needs. In the development of these products primary focus is on absolute safety to the patient achieved through ergonomic design and ease of use.

Please contact us for any additional information on any of the evomed products.

evopneu Insufflation Needles

The EVOPNEU Needles are used for traditional CO2 insufflation techniques for the simple, safe and effective creation of a pneumoperitoneum prior to abdominal endoscopy.

We offer a choice of a 2.1mm diameter needle in 120mm or 150mm lengths.



HIGH PERFORMANCE SHARP BEVELLED NEEDLE

150cm LENGTHS

RED SAFETY INDICATOR

SMOOTH STYLET RECOIL SPRING

FEATURES	ADVANTAGES	BENEFITS
150cm lengths	Thicker abdominal walls	Provides a solution for bariatric procedures
High performance sharp bevelled needle	Perforates and glides through tissue easily	Allows safe controlled insertion
Red safety indicator	Blunt safety stylet is engaged and needle rendered safe	Safety is assured and internal visceral damage avoided
Smooth stylet recoil spring	Fast safety reaction time	Minimises risk for unintentional organ damage

PRODUCT REFERENCE	DIAMETER (mm)	LENGTH (mm)
VN-120	2.1 mm	120 mm
VN-150	2.1 mm	150 mm



EVOMED Group Sarl

pH8, pépinière d'entreprises innovantes
8 place de l'Hôpital
67000 STRASBOURG, France
Phone: + 33 (0) 3 88 23 70 40
Fax: + 33 (0) 3 88 23 70 00
Web: www.evomedgroup.com



Please contact us for any additional information on any of the evomed products.